## **MODULAR HEAD-TO-TOE SOLUTION**

The PRO MAX is available with a complete selection of superior 3DEEP® RF fractional and non-fractional technologies treatment handpieces









### Intensif MX Microneedling 2x faster results

Targeted to treat acne scars, deep wrinkles, stretch marks and other conditions requiring fractional deep dermal heating.



### MiniShaper MX

Focused RF treatment for the jaw-line, submental, arms and



#### FSR MAX

Specially designed for fractional skin resurfacing and treatment of skin roughness, hyperpigmentation, wrinkles and scars on both facial and body areas.



### Shaper MX

For superior body tightening, contouring and shaping treatments



### iFine MAX

Specially designed for the treatment of wrinkles and lax skin on delicate skin tissue around the eyes and mouth



#### Contour MX

Specially designed with innovative RF + Vacuum technology to provide the ultimate in body contouring and cellulite reduction treatments.



#### Small MAX

Specially designed for skin laxity treatment of smaller areas such as the cheeks, submental, décolleté and neck.



## INFINITE POSSIBILITIES 4 TECHNOLOGIES IN ONE PLATFORM

The EndyMed PRO MAX uniquely bundles fractional and non-fractional technologies in one versatile system.



#### **Focused RF**

Energy is delivered deeper into the skin for unparalleled deep volumetric heating and clinically proven results.



The only RF technology utilizing 6 RF generators that work simultaneously, enabling complete phase control.



#### Fractional Skin Resurfacing

The only Applicator on the market with both Fractional Ablation & 3DEEP RF Heating.

Smooth scanning technology for absolute ablation control.

Ablation only effects the epidermal layer yet with superior volumetric heat penetration.

### Vacuum-Assisted Technology

Modulated pulsed vacuum technology uniquely enables deeper heat penetration for superior results - without the need for a cooling system.



## **ENDYMED**

RF Microneedling

Experience the next generation of

RF microneedle technology, now

comfort and optimal results.

with double the speed for enhanced



+972-4-630-9101 | info@endymed.com | www.endymed.com 4400 Route 9 South Freehold | Suite #1000, NJ 07728 USA

# ENDYMED PRO MAX

Superior RF Multi-Treatment Platform



## ENDYMED PRO MAX

### Superior RF Multi-Treatment Platform

The NEW 100W EndyMed PRO MAX, powered by superior, clinically validated 3DEEP® technology, offers aesthetic PRO MAXfessionals a modular head-to-toe solution that grows with your practice. Our all-in-one treatment platform is flexible, upgradeable and generates up to 7x ROI! With 3DEEP® technology, safe and suitable for all skin types and tones, you can broaden your treatment offerings and increase revenue potential with shorter appointments and effective, versatile treatments.





### **Faster** x2 microneedling

The Intensif RF Microneedling, now **x2** as fast, speeds up the treatment process, enabling you to efficiently carry out treatments and have full control over patient outcomes.



### power levels

Take advantage of power levels up to **100W** to achieve broader coverage and reduce treatment times.



### Thermal camera

The Thermal camera monitors heat levels during treatment, achieving uniform heat distribution and optimal temperatures quicker for enhanced results.



#### Versatile protocol choices


Choose from a variety of **new** treatment area protocols.





### 15-inch Intuitive GUI & new tutorials

Our intuitive, user-friendly GUI makes it easy to browse more settings and enhance your knowledge with advanced protocols and video tutorials.

### ENDYMED 3DEEP® VS. RF COMPETITORS









Monopolar RF - energy flows uncontrolled in the tissue, requiring intense cooling.

Bipolar RF - relatively superficial and less efficient.



3DEEP RF - focused and contained deep energy flow with no need for cooling.

### **CONFIDENCE IS 3DEEP®**

After 6 treatments

After 3 sessions



Dr. Stefano Toschi, Intensif & FSR



Before

After 10 sessions

After 3 sessions

Dr. Stefano Toschi, Intensif & FSR







Dr. Atul Kathed, Intensif & FSR After 3 sessions



Fresh Images Skin & Beauty Clinic, Dublin, Ireland, iFine & small After 1 treatment

### THE ENDYMED **DIFFERENCE**

#### A premium brand:

ENDYMED Medical is a worldwiderecognized brand with an extensive global presence.

### Our partnership advantage: Exceptional service, flexible upgrades and comprehensive training.

Extensive marketing support: Collateral to help your practice

grow with a high return on investment from the very start.



### **OVER 15 TREATMENT PROTOCOLS**

► Circumferential reduction

► Traumatic scar reduction

► Acne and acne scar reduction

- ► Fractional resurfacing
  - ▶ Cellulite reduction
- ► Skin tightening
- ▶ Wrinkle reduction
- ▶ Body contouring
- ► Texture improvement ► Stretch mark reduction
- ▶ Periorbital tightening
  ▶ Hyperhidrosis

► Skin rejuvenation

- ▶ Pore minimization
- ► Enhancing skin tone
- Drug delivery
- ► Facial sculpting
  - ▶ Alopecia

▶ Melasma

- ► Topical combination protocols
- \* According to treatment protocols conducted by independent physicians

As Featured In:









Los Angeles Times

